

# Inserting your earmould

Department of Audiology & Audiological Medicine  
Lanesborough Wing, St George's Hospital

This leaflet is for people who use hearing aids with an earmould. It provides step-by-step instructions for putting the earmould into your ear correctly.

## How to ensure the mould is correctly inserted into your ear

1



Hold the earmould at the back with your finger and thumb. Pull it back past your ear

2



Place the part going down the ear into position.

3



Ensure you put in place the top part of the earmould, underneath the crease at the top of your ear.

4



If you need to, use your other hand to pull your lobe out and ensure the whole of the mould is firmly in place.

5



Correctly fitted

Incorrect



Here the top part of the hearing aid is not correctly inserted under the crease of the ear. This can cause discomfort and whistling from the hearing aid.

**If you have any questions about the information in this leaflet, you can contact:**

Department of Audiology & Audiological Medicine  
Ground floor, Lanesborough Wing, St George's Hospital, Blackshaw Road, Tooting, SW17 0QT

Telephone: 020 8725 1151 (enquiries)  
e-mail: [stgh-tr.audiology@nhs.net](mailto:stgh-tr.audiology@nhs.net)  
minicom: 020 8725 3019  
text: 07552 612 845



August 2012  
Review August 2015